

VAHTS DNA Adapters Set 3 - Set 6 for Illumina

N805/N806/N807/N808

Version 24.1



Product Description

VAHTS DNA Adapters Set 3 - Set 6 for Illumina are four special kits designed for DNA library preparation for Illumina platforms. The Set 3 (Vazyme #N805) contains 24 kinds of indexed Adapters (DNA Adapter from 96-01 to 96-24). The Set 4 (Vazyme #N806) contains 24 kinds of indexed Adapters (DNA Adapter from 96-25 to 96-48); The Set 5 (Vazyme #N807) contains 24 kinds of indexed Adapters (DNA Adapter from 96-49 to 96-72); The Set 6 (Vazyme #N808) contains 24 kinds of indexed Adapters (DNA Adapter from 96-73 to 96-96). All the reagents provided in the kit have undergone rigorous quality control and functional testing to ensure the optimal stability and repeatability of library preparation.

Components

Components	N805	N806	N807	N808
From DNA Adapter 96-01 to DNA Adapter 96-24 (24 kinds)	20 µl each	----	----	----
From DNA Adapter 96-25 to DNA Adapter 96-48 (24 kinds)	----	20 µl each	----	----
From DNA Adapter 96-49 to DNA Adapter 96-72 (24 kinds)	----	----	20 µl each	----
From DNA Adapter 96-73 to DNA Adapter 96-96 (24 kinds)	----	----	----	20 µl each

Storage

Store at -30 ~ -15°C and ship at ≤0°C.

Applications

VAHTS DNA Adapters Set 3 - Set 6 for Illumina is a specialized kit for Vazyme ND series library preparation kits and suitable for multiple single-index DNA libraries for Illumina sequencers.

Notes

1. The amount of DNA Adapter used depends on the amount of input DNA.
2. Adjust the concentration of DNA Adapter to optimize for different input amount. The actual concentration of DNA Adapter refer to the instruction.

For Research Use Only. Not for use in diagnostic procedures.